

## Click Medical Hydrofluoric Acid Antidote Kit



Product Code

CM0209

Available in Sizes

**ONE SIZE**

Available in Colours



Red  
CM0209

- The National Poisons Information Service (NPIS) advise that for responsive patients, immediate decontamination is the priority for individuals exposed to hydrofluoric acid.
- Even relatively low strength hydrofluoric acid is toxic and has the potential to cause serious tissue damage and death when contaminated individuals are not treated swiftly.
  - C-Gel is a high performance calcium gluconate gel.
  - C-Gel combines with and neutralises the fluoride ion present in HF acid. This helps to prevent pain and tissue damage caused by HF acid whilst a patient is in transit to hospital.
  - C-Gel is an easy to use cream, which can be self-administered or applied by a suitable first aider, available as a 40g tube

Some examples where Hydrofluoric acids may be encountered are industrial chemistry, mining, refining, glass finishing, silicon chip manufacturing, and sanitation.

Kit contains

- 3 x C- Gel 40g
- 6 x Nitrile gloves in a carton

