

Product Code CM7003

DATASHEET

SAFETY MAT

SAFETY CHAIR

The wall mounted Safety Mat Evacuation Sled allows the smooth transfer of disabled patients over any surface even concrete or stairs. The sled can be stored in the wall mounted bag provided and easily deployed when needed in an emergency. Supplied complete with wall mounted rapid deployment pouch

- Dual inner foam provides protection and comfort for patient
- Maximum carrying weight 300kg
- 2 Safety belts and foot pocket to hold the patient safely and securely in place
- Tough slippery coated underside to assist transfer over floor surfaces and stairs
- Complete with storage pouch which can be wall mounted for rapid deployment (fixings not included)
- Dimensions (LxWxD) 1900 x 600 x 40mm

